

Module Carts

Flexibility in design & specification

MedStor[®]
The heart of your storage solutions



Healthcare storage, logistics
and product handling



MODULE CARTS FROM MEDSTOR

Bedside, Resus, Isolation, Epidural, Anaesthesia...



ORGANISATION - IDENTIFICATION - TRANSPORTATION

Multi-functional mobile storage units designed for your working environment.

CONTENTS

4-5	How to specify
6-7	Resus Cart - A mobile unit for emergency situations
8-9	Sepsis Cart - Supporting the sepsis six
10-11	Isolation Cart - Protecting against the spread of infection
12-13	Endoscope Cart - Safe transportation of surgical instruments
14-15	Trauma Cart - Essential care in life-threatening situations
16-17	General Cart - Bringing workspace and storage to any area
18-19	Anaesthesia Cart - Daily use in a busy environment
20-21	Paediatric Cart - Colourful way to deliver essential equipment
22-23	Module cart specifications
24-25	Size options - Drawer front Module Carts
26-27	Size options - Mobile Roller Shutter Module Carts
28-29	Trays and dividers
30-31	Cart security systems
32-33	Customise with accessories
34-35	Colour options & VR 360°



CONTINUOUS WORKFLOW

Module carts integrate with our healthcare storage range by utilising the Medstor tray and divider system. Medstor module carts offer a choice of:

- Height options
- Drawer configurations
- Door & drawer options
- Locking systems
- Expansive range of accessories
- Optional colourways & decals

Designed and manufactured in the UK by Medstor, giving us total system control with product continuity.



MODULE CARTS

Delivering care exactly where it is needed.

Medstor's mobile Module Carts offer flexibility in design, specification and style. They can be fully customised for purpose, so are suitable for every area from bedside to anaesthesia, built precisely for the environment in which they will be used. Complete with a choice of locks to guarantee secure storage.

Bright, robust and easy to manoeuvre, the Module Carts provide safe and efficient storage space with hygienic work surfaces and optional accessories.


The carts are compatible with Medstor's unique in-built organisation and product identification tray and liner system, which maximises storage capacity and ensures contents are easy to identify and access.

Module Carts take essential items exactly where they are needed, quickly and easily.



"A stylish, flexible way to store and transport essential items, designed around your specific needs"

Just 3.8kg/m² at 3mm




LIGHT-WEIGHT, STRONG & HYGIENIC MODULE CART PANELS

Equal rigidity to 2.4mm Aluminium (6.5kg/m²) or 1.7mm steel (13.3kg/m²)

KEY FEATURES:

Multi-functional mobile storage units designed for your working environment.

- 
CONFIGURE
- 
CONTROL
- 
CLEAN
- 
WORKFLOW
- 
GUARANTEE
- 
ACCESSORIES
- 
LIGHT-WEIGHT



BUILD TO YOUR EXACT
CONFIGURATION



EASY TO MANOEUVRE, LIGHT-
WEIGHT DESIGN



EASY TO CLEAN WORK
SURFACES



CONTINUOUS WORKFLOW
STORAGE



5 YEAR PRODUCT
GUARANTEE



WIDE RANGE OF ACCESSORIES
See pages 32-33



1025mm high (Worktop height 975mm)



Additional configurations available.

Resus Cart

A MOBILE UNIT FOR EMERGENCY SITUATIONS – TAKING LIFE-SAVING EQUIPMENT TO THE POINT OF CARE

When time is critical, having everything in one place streamlines the delivery of essential care. Medstor's mobile Resus Cart provides a platform for the equipment and consumables needed during in-hospital resuscitation, custom designed to support outstanding patient care when space is limited, but treatment vital.

"1025mm under worksurface and 1160mm bedside height options"

The Resus Cart can be built to fit available space exactly, so is always easy to store and access when needed, fully stocked and ready for use thanks to Medstor's internal tray system.

Users can specify a wide range of accessories for holding larger items such as defibrillators, oxygen cylinders and IV Poles. Accessories are securely fitted to an integral rail, so essential items are all in one place and in easy reach when needed.

CART 1025 (H) x 730 (W) x 530 (D) mm

Code	Drawer Configuration		
MC6525-2013	1 x 100mm trays	3 x 200mm trays	
MC6525-2032	3 x 100mm trays	2 x 200mm trays	
MC6525-2122	1 x 50mm tray	2 x 100mm trays	2 x 200mm trays
MC6525-2231	2 x 50mm trays	3 x 100mm trays	1 x 200mm tray
MC6525-2050	5 x 100mm trays		


CART 1160 (H) x 730 (W) x 530 (D) mm

Code	Drawer Configuration		
MC6500-2060	6 x 100mm trays		
MC6500-2013	1 x 100mm tray	3 x 200mm trays	
MC6500-2032	3 x 100mm trays	2 x 200mm trays	
MC6500-2122	1 x 50mm tray	2 x 100mm trays	2 x 200mm trays
MC6500-2131	1 x 50mm trays	3 x 100mm trays	1 x 200mm tray
MC6500-2060	6 x 100mm trays		

Please contact Medstor for additional information and personalised specification.



DRAWER FRONT MODULE CART STANDARD HEIGHT OPTIONS



See page 34 for stock colour options

RECOMMENDED ACCESSORIES

	Type	Code
A	IV pole 2 hooks	MM44792011
B	I.V. pole sharps bin holder	MM44792014
C	Defibrillator platform	MC44792070
D	Tamper-evident seal	MC49380052
E	Locking bar	MC44792052
F	Module cart handle - 452mm	MC44792040
G	Accessory rail x 2	MC44792034
H	Universal sharps bin bracket	MM44792009
I	Wally bin 8 litres	MM44792028
J	Cylinder holder	MC44792085

See pages 32-33 for full list of accessories.

Accessories are securely fitted to accessory rail located on the short side of the module cart.



RECOMMENDED ACCESSORIES

	Type	Code
A	Universal sharps bin bracket	MM44792009
B	Accessory rail x 2	MC44792034
C	Cylinder holder	MM44792085
D	Push through work surface	MC44792062
E	Module cart handle	MC44792040
F	A4 file holder	MM44794028

See pages 32-33 for full list of accessories.
Accessories are securely fitted to accessory rail located on the short side of the module cart.

Trays available in clear polycarbonate and grey ABS
see pages 28-29 for details.





1160mm high (Worktop height 1110mm)



Drawer labels clearly identified for treatment

Sepsis Cart

THE MEDSTOR CART SUPPORTS "THE SEPSIS SIX"

Ensuring that all the goods for the six stages of procedures are held together and always immediately available. Carts can be locked, see pages 30-31 for lock options.

"Supports national sepsis guidelines"

OXYGEN
BLOOD CULTURES
INTRAVENOUS ANTIBIOTICS

INTRAVENOUS FLUIDS
HAEMOGLOBIN & LACTATES
URINE MEASUREMENT



See page 34 for
stock colour options

CART 1025 (H) x 730 (W) x 530 (D) mm

Code	Drawer Configuration		
MC6525-2241	2 x 50mm trays	4 x 100mm trays	
MC6525-2330	3 x 50mm trays	3 x 100mm trays	
MC6525-2411	4 x 50mm trays	1 x 100mm trays	1 x 200mm tray
MC6525-2501	5 x 50mm trays	1 x 200mm tray	

CART 1160 (H) x 730 (W) x 530 (D) mm

Code	Drawer Configuration		
MC6500-2060	6 x 100mm trays		
MC6500-2231	2 x 50mm trays	3 x 100mm tray	1 x 200mm tray
MC6500-2312	3 x 50mm trays	1 x 100mm trays	2 x 200mm tray
MC6500-2321	3 x 50mm trays	2 x 100mm trays	1 x 200mm tray

Please contact Medstor for additional information and personalised specification.

ROOM
EMERGENCY CARE
INFORMATION RECEPTION



Isolation Cart

DESIGNED TO DELIVER CARE WHILE PROTECTING AGAINST THE SPREAD OF INFECTION

Patients with certain conditions may be nursed in isolation to prevent them spreading or contracting infection. Medstor's mobile Isolation Carts are constructed to users' specifications from products specially treated to reduce the risk of pathogen transmission and withstand rigorous cleaning.

"Having everything in one place means staff do not have to risk breaching barrier nursing protocols"

The Isolation Cart can be stocked with all necessary consumables then easily wheeled into the isolation area when needed. Medstor's tray and liner system ensures internal space is maximised and all contents are easy to see and access.

In addition to carrying consumables inside the Isolation Cart, users can also choose from a wide range of accessories for mounting larger items on the cart itself, such as glove boxes, caddie waste bins and IV Poles. Accessories are securely fitted to an external rail, so essential items are all in one place and within easy reach when needed.

CART 1025 (H) x 730 (W) x 530 (D) mm

Code	Drawer Configuration		
MC6525-2131	1 x 50mm trays	3 x 100mm trays	1 x 200mm tray
Additional configurations available.			

CART 1160 (H) x 730 (W) x 530 (D) mm

Code	Drawer Configuration		
MC6500-2231	2 x 50mm trays	3 x 100mm trays	1 x 200mm tray
Additional configurations available.			

RECOMMENDED ACCESSORIES

	Type	Code
A	IV pole 2 hooks	MM44792011
B	Caddie waste bin 7 litres	MM44792017
C	Module cart handle - 452mm	MC44792040
D	Glove box holder	MM44794003
E	Accessory rail	MC44792034
F	Apron dispenser	MM44794009

See pages 32-33 for full list of accessories.

Accessories are securely fitted to accessory rail located to the short side of the module cart.

1200mm height - MC7502 (see page 26 for tray options)

1450mm height - MC7504 (see page 27 for tray options)

1710mm height - MC7507 (see page 27 for tray options)



See page 34 for stock colour options





Endoscope Cart

FOR THE SAFE STORAGE AND TRANSPORT OF DELICATE SURGICAL INSTRUMENTS

Medstor's Endoscopy Cart is an ideal way to transport scopes and all other equipment associated with endoscopic procedures safely and efficiently, bringing vital equipment directly to the point of care and streamlining procedures.

"All designed to combine strength with enhanced infection control, essential in the surgical environment"

CART 1160 (H) x 730 (W) x 530 (D) mm

Code	Drawer Configuration		
MC6500-2231	3 x 100mm trays	2 x 200mm tray	
MC6500-2032	2 x 50mm trays	3 x 100mm trays	1 x 200mm tray

Additional configurations available.



See page 34 for stock colour options






RECOMMENDED ACCESSORIES

Type	Code
A Module cart handle - 452mm	MC44792040
B Work surface extension drop down	MM44793007
C Accessory rail	MC44792034
D IV pole 2 hooks	MM44792011
E Endoscope holder clear	MC44792056
F Endoscope holder closed (Req x 2 rails)	MC44792055
G Overbridge unit 9 bins	MC44792048

See pages 32-33 for full list of accessories.

Accessories are securely fitted to accessory rail located to the short side of the module cart.



Users can specify construction materials from a wide range of options. In addition to consumables carried inside the Endoscopy Cart, users can also choose from a wide range of accessories for mounting larger items on the cart itself, including endoscope carriers, an overbridge unit and IV poles.



CART 1025 (H) x 730 (W) x 530 (D) mm

Code	Drawer Configuration		
MC6525-2240	2 x 50mm trays	4 x 100mm trays	
MC6525-2131	1 x 50mm tray	3 x 100mm trays	1 x 200mm trays
MC6525-2050	5 x 100mm trays		
MC6525-2022	2 x 100mm trays	2 x 200mm trays	
MC6525-2003	3 x 200mm tray		

Additional configurations available.

CART 1160 (H) x 730 (W) x 530 (D) mm

Code	Drawer Configuration		
MC6500-2060	6 x 100mm trays		
MC6500-2013	1 x 100mm tray	3 x 200mm trays	
MC6500-2032	3 x 100mm trays	2 x 200mm trays	
MC6500-2122	1 x 50mm tray	2 x 100mm trays	2 x 200mm trays
MC6500-2231	2 x 50mm trays	3 x 100mm trays	1 x 200mm tray

Additional configurations available.

Trauma Cart

ENABLING THE DELIVERY OF ESSENTIAL CARE IN LIFE-THREATENING SITUATIONS


When a patient arrives at hospital following a major trauma, they need immediate access to the best diagnostic and treatment facilities, essential in life-threatening situations. Clinical staff can quickly position Medstor's Trauma Cart exactly where it is needed, with all necessary equipment in one place to avoid delays.

"Compact, manoeuvrable, robust and ideal for use in busy or limited spaces"

Trauma Cart makes up Medstor's internal tray and liner system, so can be stocked with consumables that are easily to locate and access, even in the most pressurised situation. It is also easily cleaned.

In addition to carrying consumables inside the Trauma Cart, users can also specify a wide range of accessories for mounting larger items on the cart itself, such as glove boxes, caddie waste bins and IV Poles. Accessories are securely fitted to an external rail, so essential items are all in one place and within easy reach when needed.

DRAWER FRONT MODULE CART STANDARD HEIGHT OPTIONS




See page 34 for
stock colour options

RECOMMENDED ACCESSORIES

Type	Code
A Module cart handle - 452mm	MC44792040
B Universal sharps bin bracket	MM44792009
C Glove box holder	MM44794003
D IV pole 2 hooks	MC44792011
E Overbridge unit 4 bins	MC44792050
F Accessory rail (Req x 2 rails)	MC44792034
G Bougie / Airway holder	MM44794008
H Apron dispenser	MM44794009

See pages 32-33 for full list of accessories.
Accessories are securely fitted to accessory rail
located to the short side of the module cart.





General Cart

FLEXIBLE AND EASY TO MANOEUVRE – BRINGING WORKSPACE AND STORAGE TO ANY AREA

In any hospital, there are a number of items that need to be kept together for ease of access, or to support patient care. Medstor's General Cart offers ultimate flexibility and can be used across the hospital in any department or ward to provide safe, secure storage of items and equipment and an extra work surface.

The mobile cart can also be dedicated to a particular patient, following the patient's journey and making the updating of vital information much easier and less susceptible to errors.

Medstor's General Carts are available in a wide range of bright colours; They can also be printed with any graphics that might complement surroundings, such as cartoon characters for on a children's ward.

In addition to carrying consumables inside the General Cart, users can also specify a wide range of accessories for mounting larger items on the cart itself, such as glove box holders, caddie waste bins and even a work surface extension. Accessories are securely fitted to an external rail, so essential items are all in one place and within easy reach.

"Colour-coded for ease of identification in busy hospital departments"

CART 1025 (H) x 730 (W) x 530 (D) mm

Code	Drawer Configuration		
MC6525-2013	1 x 100mm tray	3 x 200mm trays	
MC6525-2032	3 x 100mm trays	2 x 200mm trays	
MC6525-2122	1 x 50mm tray	2 x 100mm trays	2 x 200mm trays
MC6525-2231	2 x 50mm trays	3 x 100mm trays	1 x 200mm tray

CART 1160 (H) x 730 (W) x 530 (D) mm

Code	Drawer Configuration		
MC6500-2060	6 x 100mm trays		
MC6500-2013	1 x 100mm tray	3 x 200mm trays	
MC6500-2032	3 x 100mm trays	2 x 200mm trays	
MC6500-2122	1 x 50mm tray	2 x 100mm trays	2 x 200mm trays
MC6500-2231	2 x 50mm trays	3 x 100mm trays	1 x 200mm tray

Additional configurations available.



See page 34 for
stock colour options




RECOMMENDED ACCESSORIES

Type	Code
A Module cart handle - 452mm	MC44792040
B Universal sharps bin bracket	MM44792009
C Glove box holder	MM44794003
D Apron dispenser	MM44794009
E Work surface extension drop down	MC44793007
F Accessory rail (Req x 2 rails)	MC44792034
G Caddie waste bin 7 litres	MC44792017

See pages 32-33 for full list of accessories.

Accessories are securely fitted to accessory rail located to the short side of the module cart.



1200mm height - MC7502 (see page 26 for tray options)

1450mm height - MC7504 (see page 27 for tray options)

1710mm height - MC7507 (see page 27 for tray options)



Anaesthesia Cart


BUILT TO WITHSTAND THE RIGOURS OF DAILY USE IN A BUSY ENVIRONMENT

Medstor's Anaesthesia Cart is designed to store the tools needed during surgical procedures that require administration of anaesthetics. The Anaesthesia Cart helps to keep all the necessary equipment safely stored until needed, then easily at hand so they can be accessed quickly when required.

"It is vital that anaesthesia tools are well organised and maintained"

The Anaesthesia Cart is compact and easy to manoeuvre, and can be built to users' precise specifications to suit the available space and the area in which it will be used. It is also easy to clean and offers enhanced infection control.

In addition to carrying consumables inside the Anaesthesia Cart, users can also specify a wide range of accessories for mounting larger items on the cart itself, such as an overbridge unit, sharps bins and IV Poles. Accessories are securely fitted to an external rail, so essential items are all in one place and within easy reach when needed.



CART 1160 (H) x 730 (W) x 530 (D) mm

Code	Drawer Configuration		
MC6500-2060	6 x 100mm trays		
MC6500-2013	1 x 100mm tray	3 x 200mm trays	
MC6500-2032	3 x 100mm trays	2 x 200mm trays	
MC6500-2122	1 x 50mm tray	2 x 100mm trays	2 x 200mm trays
MC6500-2231	2 x 50mm trays	3 x 100mm trays	1 x 200mm tray

Additional configurations available.



See page 34 for
stock colour options

RECOMMENDED ACCESSORIES

Type	Code
A Laptop / Notebook holder	MC44792060
B Sharps bin holder	MC44792009
C Push through surface extension	MC44792062
D Wally bin 8 litres	MC44792028
E Overbridge unit 9 bins	MM44792048
F Accessory rail	MC44792034
G I.V. pole sharps bin holder	MM44792014

See pages 32-33 for full list of accessories.
Accessories are securely fitted to accessory rail
located to the short side of the module cart.



Bespoke printed graphics fit to the cart panels and brighten any paediatric department.



Broselow paediatric cart shown with overbridge unit and work surface extension.

Paediatric Carts

A COLOURFUL WAY TO DELIVER ESSENTIAL EQUIPMENT

In busy, lively paediatric wards, where health and safety may be a particular issue, the Medstor mobile Paediatric Cart is a colourful and stable way to transport essential equipment and consumables, as well as providing a useful additional work surface. Its robust construction and rounded corners make it ideally suited to a children's ward.

"Bespoke configurations and designs for your healthcare environment"

The Medstor Broselow cart is designed to be used in conjunction with The Broselow Tape (Broselow paediatric emergency tape) which relates a child's height as measured by the tape to his/her weight to provide medical instructions, including medication doses, the size of equipment that should be used and the levels of shock voltage when using a defibrillator.

All the Medstor drawers are clearly marked with the corresponding weights. Ensuring all the correct products are easily available for immediate use. Carts can be locked using tamper evident locks or dual keypad/proximity access.

Please note non-stock sizes and finishes are subject to extended lead times.




See page 34 for stock colour options



Bespoke printed graphics fit to the cart panels and brighten any paediatric department.

MODULE CART SPECIFICATIONS

- Built to your exact configuration
- Industry leading security options (See page 30)
- Full accessory options (See pages 32-33)




- | | | | |
|----|---|----|--|
| 1 | Secure fitting over bridge unit with 9,5 and 4 bin configurations. | 20 | Sharps bin holder. |
| 2 | Fully adjustable tablet holder fixes securely to cart profile. | 21 | Caddy waste bracket and bin - 7 litre. |
| 3 | Hook on glove box holder in 5mm acrylic or stainless steel. | 22 | Anti-ligature flush pull handles. |
| 4 | Seamless acrylic capped ABS top (max load 45kgs). | 23 | Reconfigure with tray dividers, designed to protect their contents. |
| 5 | Adjustable acrylic defibrillator platform holder. | 24 | Composite light weight aluminium panels. |
| 6 | Stainless steel adjustable IV pole with double loop hook. | 25 | Polycarbonate tray -40°C to plus 120°C temperature range. |
| 7 | Extend work top area with pull out and fold way surfaces. | 26 | Telescopic runners allow full withdrawal of cradle and tray. |
| 8 | I.V. sharps bracket. | 27 | Slotted dividers available in short or long sizes to fit requirements. |
| 9 | Cart universal fixing brackets and multi purpose hooks. | 28 | 4 non-marking polyolefin polymer alloy bumpers. |
| 10 | Digital locking systems include - Keypad, proximity and key options. | 29 | Rail fitted Bougie/airway holder in 5mm acrylic silver grey. |
| 11 | Aluminium profiles for strength/stability of cart panels and accessories. | 30 | Product identification cards with optional RFID dual bin. |
| 12 | Smooth oval handle rail fits to the "T" slotted profile bars. | 31 | Fully enclosed protective aluminium base panel. |
| 13 | Apron roll dispenser - hooks onto accessory rail bar. | 32 | Grey ABS tray -20°C to plus 60°C temperature range. |
| 14 | Trays are protected in a mild steel cradle. | 33 | Swivel and swivel with lock castors, high grade materials. |
| 15 | Locking bar accessory rail with tamper-evident seal. | 34 | Cylinder holder fixing clamp. |
| 16 | E23 ultrasonic welded slat - roller shutter front in grey. | 35 | Prevent static build-up - anti-static castor fitted as standard. |
| 17 | Wally bin 8 & 20 litres. | 36 | Fully welded 30mm x 30mm x 1.5mm box section. |
| 18 | Adjustable height accessory rail - fits securely to profile bars. | 37 | Fitted with 125mm Tente castors manufactured to EN 12530. |
| 19 | Laminated MDF drawer front with stainless steel handles. | | |

LIGHT-WEIGHT, STRONG & HYGIENIC MODULE CART PANELS



Just 3.8kg/m² at 3mm



See page 34 for
stock colour options

Equal rigidity to 2.4mm Aluminium (6.5kg/m²) or 1.7mm steel (13.3kg/m²)

Pencil Hardness	2H
Toughness of Coating	3T
Temperature Resistance	-50°C to +80°C
Impact Strength (kg/CM2)	50
Boiling Water Resistance	Boiling for 2hrs without change
Acid Resistance	Surface immersed in 2% HC1 for 24hrs without change
Alkali Resistance	Surface immersed in NaHO engine oil for 24hrs without change
Oil Resistance	Surface immersed in 20# engine oil for 24hrs without change
Solvent Resistance	Cleaned 100 times with Dimethylbenzene without change
Cleaning Resistance	>1000 times without change
Peel Strength	>5 Newton/mm

PANEL THICKNESS (MM)	3
Sound Absorption	0.05
Sound Attenuation	24
Water Absorption % by Volume	0.01
Thermal Performance R Values	0.0080
Core Composition	LDPE
Fire Performance	Class 1 BS476 PT7 Class 0 BS476 PT6
Tensile Strength of Aluminium Layers	149MPa
PANEL THICKNESS (MM)	-0 + 0.2
Width (mm)	±2
Length (mm)	±3
Diagonal	±5
Thermal Expansion	2.4mm/m @100°C Temp Difference
Aluminium Thickness	As specified in EN485-4

Profiles:

Lightweight anodised aluminium profiles 48mm x 22.5mm with a general wall thickness of 1.5mm are designed to secure panels from 6mm to 18mm thickness enabling the cart to be configured using: light weight composite aluminium / polymer panels.

Base frame:

Fully welded mild steel 30mm x 30mm x 1.5mm rectangular box section, electrostatic coating, with fully enclosed protective base panel.

Drawer sides:

Drawers are complete with a pair of Accuride steel telescopic runners fitted to formed and coated mild steel tray support cradles, which carry the system baskets/trays. The telescopic runners hold up to a maximum weight loading of 45kgs. Telescopic runners allow full withdrawal of the basket/trays to enable filling to the rear.

Buffers:

The Medstor Module Cart is equipped with four non marking Polyolefin polymer alloy bumpers clipped and screwed onto the fully welded base structure.

Castors:

All mobile units are fitted with 125mm Tente castors manufactured to EN 12530. Castors are chosen for roll performance, movement noise, attrition and corrosion resistance.

Swivel and swivel with lock, housing and wheels made from high grade synthetic materials, double ball bearing swivel head, bolt hole. Wheel centre made of Polypropylene, Tread: TENTEprene (thermoplastic rubber), grey non-marking, precision ball bearing. Non-static including 1 x anti-static castor


WORKTOPS

Seamless, robust and recessed acrylic capped ABS top can be fitted on as standard (max load 45kgs).



Size options

SAFE DELIVERY OF ESSENTIAL ITEMS



Drawer front module carts 1025 (H) x 730 (W) x 530mm (D)

MOBILE UNIT - GETTING EQUIPMENT TO THE POINT OF CARE

Easy to manoeuvre Tente 125mm swivel, braked and anti-static castors, buffered wheels with locking brakes. Built to house Medstor trays and dividers (see page 29).

- **Lightweight and robust cart provides hygienic working surface**
- **Easy to manoeuvre**
- **Easy tray withdrawal - full integration between trays and dividers**
- **Additional drawer configurations available on request**
- **Customise with accessories (see pages 32-33)**



See page 34 for
stock colour options



CODE: MC6525-2240
2 x 50mm trays
4 x 100mm trays



CODE: MC6525-2131
1 x 50mm tray
3 x 100mm trays
1 x 200mm tray



CODE: MC6525-2022
2 x 100mm trays
2 x 200mm trays



CODE: MC6525-2050
5 x 100mm trays



CODE: MC6525-2003
3 x 200mm tray

Additional configurations available.

Drawer front module carts 1160 (H) x 730 (W) x 530mm (D)

STORAGE AND PROTECTION - IDEALLY SUITED TO THE RIGOURS OF ANY BUSY HOSPITAL

Easy to manoeuvre Tente 125mm swivel, braked and anti-static castors, buffered wheels with locking brakes. Built to house Medstor trays and dividers (see page 29).

- **Lightweight and robust cart provides hygienic working surface**
- **Organise supplies and medication cart for bedside, isolation, epidural and anaesthesia**
- **Easy tray withdrawal - full integration between trays and dividers**
- **Additional drawer configurations available on request**
- **Customise with accessories (see pages 32-33)**



See page 34 for
stock colour options



CODE: MC6500-2060
6 x 100mm trays



CODE: MC6500-2013
1 x 100mm tray
3 x 200mm trays



CODE: MC6500-2032
3 x 100mm trays
2 x 200mm trays



CODE: MC6500-2122
1 x 50mm tray
2 x 100mm trays
2 x 200mm trays



CODE: MC6500-2231
2 x 50mm trays
3 x 100mm trays
1 x 200mm tray

Additional configurations available.



Mobile Roller Shutter Module Carts 1200 (H) x 730 (W) x 530mm (D)

MODULE ROLLER SHUTTER CARTS FEATURE & DIGITAL KEY PAD LOCK

Tambour door / Roller shutter carts with integral door braking system.

- **Welded slats with easy grip handle detail**
- **Digital locking system**
- **Additional tray configurations available on request**
- **Liner options available. Medstor ergonomic liner fitted as standard**
- **Angled or horizontal tilt arresters**
- **Rolling arresters can be fitted to ease tray withdrawal**
- **Full integration between trays and dividers**



CODE: MC7500
Trolley only



CODE: MC7500-0050
5 x 100mm ABS trays



CODE: MC7500-0031
3 x 100mm ABS trays
1 x 200mm ABS tray



CODE: MC7500-0003
3 x 200mm ABS trays

ABS TRAYS SHOWN, POLYCARBONATE ALSO AVAILABLE SEE PAGE 29



EASY TO MANOEUVRE

Tente 125mm swivel, braked and antistatic castor, buffered wheels with locking brakes.



FULLY WELDED BASE FRAME

Extremely strong yet lightweight frame design with lightweight anodised aluminium uprights in a selection of colours.



CHOICE OF DRAWER DEPTHS

Built to house the Medstor trays and dividers supplied in 50,100 and 200mm heights.

Mobile Roller Shutter Module Carts 1450 (H) x 730 (W) x 530mm (D)

MODULE ROLLER SHUTTER CARTS FEATURE & DIGITAL KEY PAD LOCK




CODE: MC7504
Trolley only



CODE: MC7504-0160
1 x 50mm ABS tray
6 x 100mm ABS trays



CODE: MC7504-0032
3 x 100mm ABS trays
2 x 200mm ABS trays



CODE: MC7504-0113
1 x 50mm ABS tray
1 x 100mm ABS tray
3 x 200mm ABS trays

ABS TRAYS SHOWN, POLYCARBONATE ALSO AVAILABLE SEE PAGE 29

Mobile Roller Shutter Module Carts 1710 (H) x 730 (W) x 530mm (D)

MODULE ROLLER SHUTTER CARTS FEATURE & DIGITAL KEY PAD LOCK




CODE: MC7507
Trolley only



CODE: MC7507-0180
1 x 50mm ABS tray
8 x 100mm ABS trays



CODE: MC7507-0133
1 X 50mm ABS tray
3 x 100mm ABS trays
3 x 200mm ABS trays



CODE: MC7507-0024
2 x 100mm ABS trays
4 x 200mm ABS trays

ABS TRAYS SHOWN, POLYCARBONATE ALSO AVAILABLE SEE PAGE 29



“The Module Carts are built around the Medstor module tray and divider system”

The Medstor trays/drawers are complete with both lengthways and widthways dividers that interwork to create multiple compartments to facilitate dual bin operations.

Trays and Dividers

STORAGE AND PROTECTION - IDEALLY SUITED TO THE RIGOURS OF ANY BUSY HOSPITAL

Extremely strong, robust and smooth Medstor trays and dividers

Above all designed to protect their contents. High temperature machine washing compliant materials are used to ensure hygienic surfaces within your environment.

Product segregation

Dividers can also be scored and snapped and used in conjunction with a system 'T' clip to create even smaller configurations. System label holders are available as standard.



TRAY TRANSPORTATION

Trays can simply be transferred from Modular Cart to any Medstor liner storage system units.

Clear Polycarbonate

TEMPERATURE RANGE -40°C TO + 120°C

- Load bearing capacity 35kg
- High impact break resistant
- Sanitised at both high and low temperatures
- Conform to international module sizes
- Full integration between trays and dividers

FULLY ADJUSTABLE TRAY DIVIDERS CREATE INDIVIDUAL COMPARTMENT

- Solid base options protect from fluid leaks
- Perforated base options improve air flow
- Inbuilt cabinet drawer stops
- Seamless hygienic design
- Nestable same size containers



EURO SIZE 600 x 400

Clear polycarbonate makes identifying materials easy, with a user temperature range of -40°C to +120°C machine wash and resistant to standard detergents.

POLYCARBONATE DIVIDERS

Correspond to relevant tray heights.

W (mm)	D (mm)	H (mm)
600	400	50
600	400	100
600	400	200



OPEN SIDES
AND SOLID BASE



OPEN SIDES
AND OPEN BASE



SOLID SIDES
AND SOLID BASE



SHORT DIVIDER
400 D (mm)



LONG DIVIDER
600 D (mm)

600	400	50	-	-	M1605
600	400	100	M1601	M1651	M1611
600	400	200	M1602	-	M1612

M1645	M1665
M1641	M1661
M1642	M1662

ABS

TEMPERATURE RANGE -20°C TO + 60°C

- Load bearing capacity 35kg
- Good impact strength
- Conform to international module sizes
- Full integration between trays and dividers

FULLY ADJUSTABLE TRAY DIVIDERS CREATE INDIVIDUAL COMPARTMENT

- Anti-static tray dissipates static charge build-up
- Solid base options protect against fluid leaks
- Perforated base options improve air flow
- Inbuilt cabinet drawer stops
- Seamless hygienic design
- Nestable same size containers



EURO SIZE 600 x 400

ABS offers a tidy container solution, with a user temperature range of -20°C to +60°C resistant to standard detergents.

ABS DIVIDERS

Correspond to relevant tray heights.

W (mm)	D (mm)	H (mm)
600	400	50
600	400	100
600	400	200



OPEN SIDES
AND SOLID BASE



OPEN SIDES
AND OPEN BASE



SOLID SIDES
AND SOLID BASE



SHORT DIVIDER
400 D (mm)



LONG DIVIDER
600 D (mm)

600	400	50	-	-	M0605
600	400	100	M0601	M0651	M0611
600	400	200	M0602	-	M0612

M0645	M0665
M0641	M0661
M0642	M0662

Cart Security Systems

SECURE STORAGE IN ANY ENVIRONMENT



Medstar Module Carts can be fitted with a variety of security locking options, ensuring that dangerous, costly or sensitive items are stored safely, wherever the Module Cart is being used.

Users can choose from digital code security pad or key locking systems, a locking bar or proximity fob options.

Fully tamper evident seals and individual drawer locks are also available.

"Industry-leading Module Cart security"

Clinical staff using Medstar Module Carts can be sure that even when they are working in busy or pressurised environments, items in the cart are safe and secure, offering invaluable peace of mind and allowing them to focus on patient care.



Digital locking

Lock code: MC49380052 Card code: MC49380063



KEY LOCKING

Suited claw lock key option -
Key opens all drawers.

Key locking

Code: MC44792063



LOCKING BAR

Cost effective cart security
used with a security seal.

Locking bar

Code: MC44792061



INDIVIDUAL DRAWER

Tamper security seal and identification
for outer cart.

Tamper evident seal

Code: MC49380052



WE UNDERSTAND THE DIFFICULTY OF TRACKING INVENTORY

Please contact Medstar for full details on our RFID digital systems inventory management solutions and transponder technology for cabinets and modular carts.





ANTI-LIGATURE HANDLES

Built into drawer fronts, available on request - stainless steel handles fitted as standard.

HANDLE OPTIONS

Stainless steel & anti-ligature

Stainless steel handles fitted as standard.

Customise with Accessories


A TAILORED APPROACH TO CART DESIGN

Medstor Module Carts are designed to maximise all available space and ensure essential items are together and within easy reach when needed. A wide range of accessories enables users to specify a bespoke design that fits their needs exactly.


Larger items or those that need to be accessed frequently like sterile gel, can be mounted on the cart itself: glove box holders, caddy waste bins and even a defibrillator platform all attach securely via an adjustable accessory rail. For maximum safety, ergonomic design ensures all items fit as flush to the cart as possible.



WORK SURFACE EXTENSION
Push through: MC44792062



WORK SURFACE EXTENSION
Drop down: MM44793007



ENDOSCOPE HOLDER
5 spider: MC44792056
Closed carrier: MC44792055
See pages 12-13



VIAL HOLDER SHELF
Clear: MC44792054



PLASTER BOWL HOLDER
6 litre stainless steel bowl
& holder: MC44792009



LAPTOP / NOTEBOOK HOLDER
Silver grey: MC44792060



TABLET HOLDER
Silver grey: MC44792080



DRAIN
LSU shelf: MC44792012



ADJUSTABLE HANDLE
Short side 452mm: MC44792040
Long side 652mm: MC44792041



TRAY EXTENSION
Code: MC44792046




CPR BOARD
5mm acrylic: MC44792013



I.V. POLE 2 HOOKS
Stainless steel: MM44792011



I.V. pole sharps bin holder
Bracket Code: MM44792014



DEFIB PLATFORM
Grey 5mm acrylic: MC44792070



ELECTRICAL ACCESSORIES
UK & Euro fittings available
Code: MC49380044



OVERBRIDGE UNIT
4 bin: MC44792050
5 bin: MC44792049
9 bin: MC44792048
With shelf: MC44792053



ACCESSORY FITTINGS

The Medstor Module Cart is complete with a "T Slot" accessory facility. A smooth flat oval handle and multiple accessory rails slide into the upright.



HANDLE TYPES

- 1** Stainless steel handles fitted as standard
- 2** Anti-ligature handles fronts available on request

ACCESSORY RAIL MOUNTED ACCESSORIES



DIN ACCESSORY RAIL

Short side 474mm: MC44792034

Drawer front long side: MC44792029

Roller shutter long side: MC44792031



CYLINDER HOLDER

Code: MM44792085



GLOVE BOX HOLDER

Grey 5mm acrylic: MM44794003



GLOVE HOLDER - TOP LOADING

Grey 5mm acrylic: MM44794007




BOUGIE / AIRWAY HOLDER

Grey 5mm acrylic: MM44794008



CADDIE WASTE BIN

7 litres: MM44792017



WALLY BIN

8 litres: MM44792028

20 litres: MM44792029



FILE HOLDER

Grey 5mm acrylic: MM44794028



SHARPS BIN HOLDER*

Black Polymer: MM44792009

*Bracket only



APRON DISPENSER

Grey 5mm acrylic: MM44794009



GEL / BOTTLE HOLDER

Grey 5mm acrylic: MM44794005

"When it comes to building a one-off design, our experienced team is on hand to help"

STOCK COLOUR PALETTE Additional colour options available on request.



COLOUR
Light Grey - GY
Egger Light Grey U708

COLOUR
Platinum White - WH
Egger Platinum White W980

COLOUR
China Red - CR
Egger China Red U321

COLOUR
Brilliant Yellow - YL
Egger Brilliant Yellow U114



COLOUR
Blue Flash - BF
Polyrey B001

COLOUR
Blue Turquoise - LB
Polyrey B119

COLOUR
Pervenche - BR
Polyrey P016

COLOUR
Mousse - EG
Polyrey M047



COLOUR
Rouge Griotte - PR
Polyrey R047

COLOUR
Cactus - CG
Polyrey C049

COLOUR
Pomme Vert - PV
Polyrey P092

COLOUR
Florida Orange - FL
Polyrey F001

Colour options

ADDITIONAL COLOURS ON REQUEST

Medstor Module Carts are available in a variety of stock and non-stock colours and finishes. If colour accuracy is important, we suggest ordering a physical swatch samples card for your Module Cart.

- **Ward colour matching service**
- **Colour coordinated panels**
- **Fitted decals**

Non-stock colours subject to extended lead times



"Material swatches are available to order from our sales team"



LOOK OUT FOR THE ICON

Are you ready to start viewing the Medstor Virtual Reality and 360° experiences?



To get started you will need a compatible VR headset or you can view our 360° experience directly from your internet website browser:

Please visit www.medstor.co.uk/vr



Require technical help?

Call us on +44 (0)161 330 5577 for help setting up this new way to experience Medstor healthcare products.



REQUIRE BIM BLOCKS?

In-house space planning team offer bespoke design service to transform your workplace environment including stock level calculations and space management.



Our BIM block collection is available to download for architects and specifiers on request



Flexibility in design, specification and **style**, multi-functional mobile units compatible with Medstor trays for **maximum organisation**.
Customise for purpose with practical accessories.

UK HEAD OFFICE & SHOWROOM
Bayley Street,
Stalybridge,
Cheshire, SK15 1QQ

TEL: +44 (0)161 330 5577
EMAIL: health@medstor.co.uk

www.medstor.co.uk

IRELAND OFFICE
Unit 120, Grange Hill
Baldoye Industrial Estate,
Dublin 13.

TEL: +353 (0)1 839 3596
EMAIL: sales@medstor.ie

www.medstor.ie