


Theatre Power Pod

Multi-functional O.R. aid

MedStor[®]
The heart of your storage solutions



SHOW EDITION
STAND NUMBER: HALL 13 A28



ILLUMINATED POWER UNITS

BRIGHTEN ANY OPERATING
THEATRE AREA WITH LED DISC LIGHTS



SAFETY INFORMATION

Theatre power pod reduces potential trip hazards
when fitted with our LED disc light.

WEIGHTED ANTI-TILT MOBILE BASE

STAINLESS STEEL TELESCOPIC
2 & 4 IV HOOKS




Code: MR1241-42-43-2

Code: MR1241-4

Code: MR1241-43-2-E

Code: MR1241-42



ALL POD POLYMER PRODUCTS MANUFACTURED FROM FLAME RETARDANT ML54 MDPE WITH ADDED ANTIMICROBIAL.

32mm diameter stainless steel IV pole with telescopic solid stainless steel rod and pig tails. Choose from 2 or 4.

UK, EURO, USA SOCKETS AVAILABLE ON REQUEST.

EASY TO MANOEUVRE

Medical graded swivel castors for high end applications, this castor increases the IV pole stability through its total lock break stem. Single anti-static caster fitted as standard.



WEIGHTED BASE AND IV POLE/DRIP STANDS TO BE TESTED TO BS3619,1979 CLASS 1 DEVICE, UNITS TO BE CE MARKED.

HYGIENIC DESIGN FUNCTIONALITY

Removable theatre pod accessories

POD ITEMS FIT SNUGGLY INTO EACH OTHER AND UNITS CAN BE MIXED AND MATCHED.

- Hygienic removable handle for easy cleaning
- Removable MDPE carry container with flame retardency and antimicrobial additive
- Illuminated power unit to be certified BS EN 60601.




ALL UNITS ARE COMPLETE WITH PUSH HANDLES



ROBUST CARRY CONTAINER 24 LITRES



UK OR EUROPEAN x5 TWIN SOCKETS WITH 3 METRE FLEX



ILLUMINATED POWER UNIT

Code	Description - Base weighted for added stability
MR1241-2	Base with 5 castors, 2 hook IV pole
MR1241-4	Base with 5 castors, 4 hook IV pole
MR1241-42	Base with 5 castors, Carry container, Pole
MR1241-42-2	Base with 5 castors, Carry container, 2 hook IV Pole
MR1241-42-4	Base with 5 castors, Carry container, 4 hook IV Pole

Code	Description - Base weighted for added stability
MR1241-43	Base with 5 castors, Illuminated Power Unit
MR1241-43-2	Base with 5 castors, Illuminated Power Unit, 2 hook IV Pole
MR1241-43-4	Base with 5 castors, Illuminated Power Unit, 4 hook IV Pole
MR1241-42-43	Base with 5 castors, Carry container, Illuminated Power Unit
MR1241-42-43-2	Base with 5 castors, Carry container, Illuminated Power Unit, 2 hook IV Pole
MR1241-42-43-4	Base with 5 castors, Carry container, Illuminated Power Unit, 4 hook IV Pole



UK HEAD OFFICE & SHOWROOM
Bayley Street,
Stalybridge,
Cheshire, SK15 1QQ

TEL: +44 (0)161 330 5577
EMAIL: health@medstor.co.uk

www.medstor.co.uk

IRELAND OFFICE
Unit 120, Grange Hill
Baldoyle Industrial Estate,
Dublin 13.

TEL: +353 (0)1 839 3596
EMAIL: sales@medstor.ie

www.medstor.ie